

Catalog Number

70000-HG6

Guard Column Cartridge, Diamond Hydride, Replacement Cartridge, 4.0mm ID x 10mm Length, 4um, 100A. In an individual black case. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [Hichrom Style Guard Column Holder](#).

Sometimes, the individual cartridge is preferred over buying a pack of 5 of the same phases for when different "Guard Column" phases are preferred for pairing them with different analytical phases for unique separations.

Use this guard column cartridge with a 3.0, 4.0 and 4.6mm ID analytical column.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$273.86

Printed On 2024/05/20