

Catalog Number

70000-SP-G1

Semi-Prep Guard Column Cartridges, Diamond Hydride,
Replacement Cartridges, 10mm ID x 10mm Long, 4um, 100A. Cogent
TYPE-C Brand. 3/PK.



Additional Info

THIS IS A CUSTOM ORDER.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,996.91

Printed On 2024/05/19