

## Catalog Number

75001-15P

**HPLC Column, HPS Silica, 5um, 4.6mm ID x 150mm Length, 120A.**

Cogent Brand. 1 EA.

**MOST POPULAR.**



### Additional Info

This Very Popular Cogent HPS Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This proprietary Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification or for new methods.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles and is often used for hydrophilic or polar compounds in Normal Phase HPLC.

*Compare its performance to Inertsil and Luna products.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,099.81**

Printed On 2024/05/20