

Catalog Number

75020-15P

HPLC Column, HPS Phenyl, 5um, 4.6mm ID x 150mm Length, 120A.
Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available making it ideal as an Alternate in USP Methods that use the L11 Classification.

With a 5µm Particle, this Column works well with higher flow rates producing lower back pressures compared to columns with smaller particles. It is often used to separate Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC.

Compare this column's performance to Luna™ or Inertsil™ columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,060.16

Printed On 2024/05/20