

Catalog Number

75025-P150

Prep HPLC Column, HPS Cyano, 5um, 120A. 21.2mm ID x 150mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent HPS Cyano™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond, and it is bonded to our 120A High Purity Silica, making this Column ideal for purification of polar compounds in Reversed Phase.

With a 5 µm Particle, this Column works well with higher flow rates and lower back pressures while the large diameter and length of this column provides high loading capacity which is great for Preparative HPLC.

Compare the performance of this Column to Inertsil or Luna brands.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$8,522.66

Printed On 2024/05/20