

## Catalog Number

75040-15P-2

**HPLC Column, HPS Amino**, 5um, 2.1mm ID x 150mm Length, 120A.  
Cogent Brand. 1 EA.



### Additional Info

This Cogent HPS Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using tri-functional bonding of a 3 carbon silane with a terminal amino group to our unique High Purity Silica. This Column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 5 µm Particle, this Column works well with higher flow rates and lower back pressures compared with smaller particles. Used for polar compounds with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or HILIC modes. Compare to Inertsil™ Amino and Luna™ Amino Columns.

*A very popular column size, this narrow-bore column provides great data but requires a pump that can precisely deliver, and a flow cell that is made for, reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,060.16**

Printed On 2024/05/20