

Catalog Number

82361-25P-4-3M

HPLC Column, NP Silica, Wide Pore, 3 μ m, 4.0mm ID x 250mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent NP Silica™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 300A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification for large molecules or for new methods.

With a 3 μ m Particle, this column is more efficient than a 5 μ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure.

Often used for large hydrophilic or polar compounds in HILIC or Normal Phase HPLC.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,036.93

Printed On 2024/05/20