

## Catalog Number

95001-01ND

**Inserts, Low Volume** for Vials, Glass. Clear, Deactivated, 200ul fill volume with a conical, precision point and 6x29mm outer dimensions. For use in 8mm, Push Top, 1ml autosampler vials. AQ Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a precision point conical tip for maximum recovery of sample. Residual volume is 2ul with most autosamplers. This product is manufactured in an ISO 9001:2015 Certified facility.

*Spring may be required, sold separately, Can be used with 8x40mm, "96 position" vials.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$160.08**

Printed On 2024/05/20