

## Catalog Number

95001-0CP-INS

**Vial & Cap Kit** includes 100 EA 1ml, Push Top, Clear Autosampler Vials and 100 EA matching, translucent, polyethylene Push Caps for Micro-Insert. BASIK Brand Easy Purchase Pack. 1/PK.



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps.

Shell Vials made with 1st Hydrolytic, Pharmaceutical grade Borosilicate glass with outer dimensions 40x8.26mm.

High purity polyethylene Push Top Caps are designed for use with 150ul Micro-Inserts [95001-INS](#) in Waters or Shimadzu autosamplers that accept 8x40 vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$40.04**

Printed On 2024/05/19