


## Catalog Number

95025-01-2S

**Caps, Headspace**, Crimp Top, with Dark Grey Butyl Rubber Septa, fitted in Aluminum, Silver Pressure Relief Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



### Additional Info

Molded, Butyl Rubber Septa is 55° Shore A (hardness) and is 3.0mm thick. Temperature resistant from -40C to 120C. Recommended for non critical analyses or when PTFE is an analyte or must be avoided.

Clear lacquered Caps have score lines that will break open at 3.0 +/- 0.5 Bar, for safety. Low cost for budget minded laboratories.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$49.09**

Printed On 2024/05/20