


Catalog Number

95025-07-1S

Caps, Headspace, Crimp Top, with White Silicone Rubber / Natural Beige PTFE Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand.



Additional Info

Silicone Rubber / PTFE with Durometer of 45 Shore A, (hardness) and 3.0mm thick. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Low cost and one of our Most Popular Headspace cap.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$86.06

Printed On 2024/05/20