


Catalog Number

95025-26-1D

Caps, Headspace, Crimp Top, with White Silicone Rubber / Blue PTFE
SPME Style Septa, fitted in Steel, Magnetic, Gold Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand.



Additional Info

Special Shape Septa with thick rubber outer ring (3mm) and thinned rubber center (1.5mm) provides a good seal with most 20mm crimp top Headspace vials. The target area on this Septa is easier for SPME needles to penetrate, minimizing broken or damaged needles..

Recommended with CTC Pal® and other autosamplers with SPME or Solid Phase Micro-Extraction injectors.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$84.59

Printed On 2024/05/20