


Catalog Number

95025-60-C

Caps, Headspace, Crimp Top, with Grey Butyl Rubber / Transparent PTFE Septa, fitted in Steel, Magnetic, Gold Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand. 1000/CS.



Additional Info

Butyl Rubber Septa that is laminated with PTFE and is 55⁰ Shore A (hardness) and is 3.0mm thick.

Temperature resistant from -40C to 120C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,494.83

Printed On 2024/05/20