

Catalog Number

95025-E-PE10A

Vials, Headspace, Crimp Top, Glass. Amber. 10ml, 20mm Serum Finish, 23x46mm outer dimensions, Beveled Top and Rounded Bottom with a Tapered Neck. Perkin Elmer Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Perkin Elmer®, Tekmar® and other Headspace autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$72.49

Printed On 2024/05/20