

Catalog Number

95025-R

Rack for Headspace Vials, Clear acrylic plastic. 25 cavities with a diameter of 24mm. MicroSolv Brand. 1 EA.



Additional Info

Often used for storage, transportation or filling, this heavy weight rack is excellent for 20mm wide Headspace Vials. The clear acrylic makes it possible to see into the bottom of the vials for problems with your samples.

Dimensions: 160mm x 160mm x 30mm.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$158.96

Printed On 2024/05/19