


## Catalog Number

9502C-0-WAV

**Vials, Crimp Top, Glass. Amber, 2ml, with Write-On Patch and fill lines.**  
An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. These wide mouth, amber vials are often used for light sensitive compounds and considered low actinic.

Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

*Use with 11mm aluminum crimp caps. These vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$44.38**

Printed On 2024/05/20