

## Catalog Number

9502C-00-S-M

**Caps, Crimp Top**, with an Open Hole. Aluminum, Silver Caps. Septa are sold separately and must be inserted manually. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand.



### Additional Info

Click [HERE](#) for Septa page for this cap. Choose Septa for 11mm Caps.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$88.04**

Printed On 2024/05/20