


## Catalog Number

9502C-10CY-M

**Caps, Snap Top**, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Yellow Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



### Additional Info

This cap is our original Snap Cap and has been on the market for 30 years. Case quantities are the most economical value.

Click [HERE](#) for smaller package quantities.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$341.34**

Printed On 2024/05/20