

Catalog Number

9502C-31-MB

Caps, Snap Top, with White Silicone Rubber / Tan PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. AQ Brand. 100/PK.



Additional Info

This is one of the most popular snap caps on the market and has been available for over 25 years. *Pre-Slit septa are normally recommended for Waters Alliance instruments.*

Compares to Waters® HPLC Certified Snap Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$70.64

Printed On 2024/05/19