


## Catalog Number

9502C-40T-R-M

**Caps, Crimp Top**, with Red Rubber / PTFE Septa, fitted in aluminum, Red Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

The Rubber is softer than Natural Red Rubber and provides better sealing after injection. Great for GC.

Click [HERE](#) for smaller pack quantities.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$198.51**

Printed On 2024/05/20