

## Catalog Number

9502C-S30-B

**Caps, Snap Top**, with 1.00mm thick, White Silicone Rubber / Blue PTFE Pre-Slit (X Shape) Septa, fitted in smooth, soft polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv



### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$87.31**

Printed On 2024/05/20