

## Catalog Number

9502R-02B-K

**Vial Racks**, for 12x32mm, 2ml autosampler vials. Polypropylene, blue, with 50 wells that are indexed. Stackable. MicroSolv brand. 5/PK.



### Additional Info

Convenient and easy management for transportation or storage of standard, 12mm wide autosampler vials that are used in HPLC, LCMS, GC or GCMS. This polypropylene rack is chemically resistant to most HPLC solvents. Racks are stackable for storing them.

Conical and flat bottom vials can be accommodated in this rack. Click [HERE](#) for single rack quantity.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$86.19**

Printed On 2024/05/19