


Catalog Number

9502S-01-B-M

Caps, Screw Top, with an Open Hole. "Knurled", polypropylene, Blue Caps. Septa are sold separately and must be inserted manually. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 1,000/CS.



Additional Info

These caps are used for our Soft-Guard Brand Caps.

Click [HERE](#) for Septa to use in these caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.69

Printed On 2024/05/19