

Catalog Number

9502S-07-LCC

Caps, Screw Top, with Clear Silicone Rubber / PTFE, ultra-pure, Pre-Slit LC-MS compatible LC Clean Septa, bonded in "knurled", polypropylene, Grey Caps. For use with 9-425mm thread, 12x32mm autosampler vials.



Additional Info

LC Clean™, LCMS compatible.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$66.19

Printed On 2024/05/20