

Catalog Number

9502S-09-SR-K

Caps, Screw Top, Solid, with a Red-Orange Silicone Rubber / TEF liner, in "knurled", closed Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand 1000/CS.



Additional Info

Silicone Rubber with TEF liner is inert lowers the cost. Liner is 53° Shore D, 0-2mm.

Solid top turns autosampler vials into storage vials. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$339.53

Printed On 2024/05/20