


## Catalog Number

9502S-09-ST

**Caps, Screw Top, Solid**, with a 0.2mm thick, PTFE only liner, in "knurled", closed Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

PTFE is a very inert material for liners and is compatible with most chromatography samples and solvents. However, PTFE can be minimally porous and does not offer the same sealing capabilities as a Rubber / PTFE liner may provide. This cap is very useful for storage of non-volatile samples.

*These solid polypropylene, blue caps turn autosampler vials into storage vials.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$43.85**

Printed On 2024/05/20