

Catalog Number

9502S-0AV-K

Vials, Screw Top, Glass. Amber, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. (10 packs of 100/PK). 1,000/CS.



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized.

Features a wide mouth opening for easy filling and ND9 threads. Lowest cost in class.

Case quantity is the most economical value.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$310.63

Printed On 2024/05/20