

Catalog Number

9502S-11-CB

Caps, Screw Top, with Natural Silicone Rubber / White PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Natural Silicone Rubber is made without adding color and has the same properties as White Silicone Rubber and is a less expensive bonded Septa and Cap.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$44.88

Printed On 2024/05/20