

## Catalog Number

9502S-35-B

**Caps, Screw Top**, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth sided, Robo Style, soft polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge" into vials. Compares to Thermo / National Scientific DP® Caps.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$39.58**

Printed On 2024/05/20