


Catalog Number

9502S-35-Y-M

Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth sided, Robo Style, soft polypropylene, Yellow Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand.



Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge" into vials. Compares to Thermo / National Scientific DP® Caps.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$328.69

Printed On 2024/05/20