

## Catalog Number

9502S-3XWCP

**Vial & Cap Kit** Includes 100 EA 2ml, Screw top, Clear Autosampler Vials with Write-On Patch and fill lines and 100 EA matching Screw Caps with Silicone Rubber / Red PTFE Pre-Slit (X Shape) Septa, fitted in the Blue



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Vials: Cat. No [9502S-WCV](#) and Caps: Cat No. [9502S-E-3XB](#).

This is our Most Popular standard borosilicate glass vial and cap Easy Purchase Pack. For our Most Popular RSA™ Vial Kits, Click [HERE](#).

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$76.10**

Printed On 2024/05/20