

Catalog Number

9502S-80C-B-M

Caps, Screw Top, with Natural Rubber / TFE Septa, fitted in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

Natural Rubber does not contain siloxanes like Silicone Rubber. TFE is similar to PTFE when comparing inertness.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$346.20

Printed On 2024/05/20