

Catalog Number

9502S-CAP1

Caps, Screw Top, with an Open Hole. "knurled", polypropylene, Blue Caps. Septa are sold separately and must be inserted manually. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard Brand.



Additional Info

Click [HERE](#) for Septa Only page.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$54.81

Printed On 2024/05/19