

Catalog Number

9502S-E-35B-M

Caps, Screw Top, with White Silicone Rubber / Natural PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

This is our most economical and lowest cost Screw Cap with Pre-Slit Septa for Autosampler Vials.

Natural PTFE is equally inert and as durable as Red PTFE, but dyes have been added to it. Normally recommended for Waters® Alliance and Acquity instruments.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$294.45

Printed On 2024/05/20