

## Catalog Number

9502S-E-35K-M

**Caps, Screw Top**, with White Silicone Rubber / Natural PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Black Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

*This is our most economical and lowest cost Screw Cap with Pre-Slit Septa for Autosampler Vials.*

Natural PTFE is equally inert and as durable as Red PTFE, but dyes have been added to it. Normally recommended for Waters® Alliance and Acquity instruments.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$297.13**

Printed On 2024/05/20