

Catalog Number

9502S-E-WAV

Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Features a wide mouth for easy filling. Low cost for budget minded laboratories.

Amber vials are often used for light sensitive compounds and considered low actinic. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Lowest cost in class for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$21.84

Printed On 2024/05/20