

## Catalog Number

9502S-PP-AMB-K

**Vials, Screw Top, Plastic. Amber, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.**



### Additional Info

Molded plastic vials are made with virgin polypropylene and have fill lines and graduations. Amber color is useful for light sensitive compounds.

Excellent choice for Ion Chromatography, CE or Bio-analysis. Minimize ionic interactions with vial walls.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$527.13**

Printed On 2024/05/20