

## Catalog Number

9502S-PP-CLR-K

**Vials, Screw Top, Plastic.** Translucent, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.



### Additional Info

Our most popular plastic autosampler vial can be used instead of glass vials for HPLC. These high volume (2ml) blow-molded, plastic vials are made with virgin polypropylene, have fill lines with graduations and a slightly concave bottom that is helpful in reducing residual volume.

Non-ionic, these PP vials can reduce or eliminate adsorption of many compounds including proteins, peptides and bases where adsorption is common with borosilicate glass.

Excellent choice for Ion Chromatography, CE, Bio-analysis, analysis of [PFAS](#) or simply replacing basic HPLC autosampler vials. Compatible with most modern autosamplers that accept 12x32mm vials.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$472.00**

Printed On 2024/06/01