

Catalog Number

9502S-WDAV

Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, Deactivated, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Deactivated / Silanized pharmaceutical grade, Borosilicate Glass, amberized. Write-On Patch and fill lines are epoxy paint for reduced smearing. Features a wide mouth opening for easy filling and ND9 threads.

Amber vials are often used for light sensitive compounds and are considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$92.49

Printed On 2024/05/20