


## Catalog Number

9502SB-06

**Storage Box** for 2ml Autosampler Vials. Natural, chemical resistant polypropylene holds 81 vials. Complete with translucent plastic lid. 1 EA.



## Additional Info

Excellent for standing position on the lab bench or during transport and ideal for space saving storage in places such as refrigerators.

- Lid prevents condensation on the caps minimizing contamination during storage.
- Usable temperature range -80C to +100C.
- Alpha numeric indexing coded on bottom of wells and along the sides of the box for clear sample identification.
- Box outer dimensions are 130mm wide by 130mm deep by 45mm high.
- Box and lid are Autoclavable.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$43.60**

Printed On 2024/05/19