

Catalog Number

9503S-E-0AV

Vials, Screw Top, Glass. Amber, 2ml, 10-425mm threads, 12x32mm
outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. . Features a wide mouth opening for easy filling and ND10 threads.

Amber vials are often used for light sensitive compounds and are considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$34.34

Printed On 2024/05/20