

Catalog Number

9504S-03ND

Inserts, Low Volume for Vials, Glass. Clear, Deactivated, 700ul fill volume with a conical, precision point and 8x40mm outer dimensions with plastic spring attached. For use in 4ml autosampler vials. AQ Brand.



Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a precision point conical tip for excellent recovery of sample. Residual volume is 8ul with most autosamplers with a 500ul useable volume. Polypropylene springs are attached to keep the inserts center aligned and reduce chances of the needle breaking the inserts during use. This product is manufactured in an ISO 9001:2015 Certified facility.

Use with 4ml, 13mm, Screw Top, 15x45mm autosampler vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$186.05

Printed On 2024/05/20