

Catalog Number

9504S-05C-BLK-M

Caps, Screw Top, with Silicone Rubber / PTFE Septa, fitted in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

Economical Cap with Septa for use in Older Waters instruments that accommodate 4ml autosampler vials.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$270.00

Printed On 2024/05/20