

Catalog Number

9504S-05S-W

Caps, Screw Top, Solid, with a Silicone Rubber / PTFE liner, in "knurled", polypropylene, White Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand 100/PK.



Additional Info

Turns 15x45mm, 4ml autosampler vials into storage bottles. Silicone Rubber / PTFE liner provides a good inert closure with good resealing properties.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$83.44

Printed On 2024/05/20