

## Catalog Number

9504S-08-SP

**Caps, Screw Top, Solid, with a Red Butyl Rubber/ Grey PTFE liner, in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand. 100/PK.**



### Additional Info

Economical Cap turns 15x45mm, 4ml autosampler vials into storage bottles. Butyl Rubber / PTFE liner provides a good inert closure with adequate re-sealing properties.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$56.60**

Printed On 2024/05/20