

## Catalog Number

9504S-0AV

**Vials, Screw Top, Glass. Amber, 4ml, 13-425mm threads, 15x45mm**  
outer dimensions. For use as an autosampler vial or storage bottle.  
MicroSolv Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized.  
Features ND13 threads.

Amber vials are often used for light sensitive compounds and considered low actinic.

*Use these vials with Waters WISP® and older Shimadzu® autosamplers.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$51.41**

Printed On 2024/05/20