

Catalog Number

9504S-10C-BK

Caps, Screw Top, with Silicone Rubber / PTFE Septa, bonded in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Use with 15x45, 4ml autosampler vials. Used with older Waters and Shimadzu autosamplers.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$82.33

Printed On 2024/05/20