

Catalog Number

9508C-10C-D

Caps, Push Top, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Natural Transparent Caps. For use with 8mm, Crimp Top autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

This cap is similar to a Snap Cap but does not need the Snap Ring on a vial to create a tight seal.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$55.89

Printed On 2024/05/20