


## Catalog Number

9509C-31CB

**Caps, Snap Top**, with Silicone Rubber / PTFE, ultra-pure, Pre-Slit (Star Shape) Septa, fitted in smooth, polyethylene, Light Blue Caps. For use with 11mm, 12x32mm, Snap Top autosampler vials. AQR Brand. 100/PK.



### Additional Info

This specially designed septa provides easy penetration by autosampler needles with improved resealing compared to other Snap Caps. The ultra-pure Silicone Rubber and PTFE have no added color producing the highest purity septa on the market with virtually no contamination or extractables.

The Star Shaped (\*) slit offers a much greater target area for an autosampler needle to penetrate the septa and offers minimal resistance. It may not reseal as well as other Pre-Slit Septa.

Recommended for LCMS/MS.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/20