

## Catalog Number

9509S-0CV-RS

**Vials, Screw Top, Glass. Clear, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. RSA Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type-33 Reduced Surface Activity Borosilicate Glass. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

Low adsorbing and will not add Na to sample solutions or compounds that cause pH changes. Click [HERE](#) for more information about RSA autosampler vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/20